



OXYGEN CONCENTRATOR

EFFICIENT AND CONVENIENT

GENE O₂ provides Oxygen Therapy for Hospital, Care Homes & Home Usage, delivers Oxygen at 93±3% Purity as per the Medical Standard.

Pressure Swing Adsorption Technology

GENE O₂

DECK MOUNT



Continuous Flow Setting

Uninterrupted Supply of Oxygen



Oxygen Concentration 93±3%

Ensuring Purity



Key Features

Gene O₂ oxygen concentrator provides supplementary oxygen for patient with COPD and in higher concentration and low flow oxygen therapy for severe hypoxemia and pulmonary edema.

GENE O₂ assures 93 ± 3% purity of Oxygen supply. Suitable for Hospitals, Nursing/ Patient/ Care Homes and Homes.

GENE O₂ uses ambient air to generate oxygen using PSA Technology.

Warranty (Excluding Sieve) - 3 Years

Product Code

Oxygen Concentrator **Gene O₂ 5 Liter/ 10 Liter**

Technical Specification

Voltage	230V+/- 10%/ 50Hz+/-1 Hz
Highest flow	10 Liter Per Min/ 5 Liter Per Min**
Power consumption	770 VA/ 350 VA **
Weight	20.5 Kg/ 17 Kg** approx
Dimensions (L*W*H}	L 370*W 300*H 540 (mm)
Oxygen output pressure	40 kPA - 68 kPA
Purity	93 ± 3%
Range of adjusting flow	0.5 to 10 Liter Per Minute/ 0.5 to 5 Liter Per Minute**
Noise	Less than <55 dBA/ <46 dBA**
Country of Origin	India
Alarms	<ul style="list-style-type: none">When purity <80%, the O₂ Low indicator will flash. This will switch off when purity is >80%.When purity < 80% for more than 3 minutes, the call service alarm will flash.In case call service indicator flashes continuously, please turn off at once, arrange spare oxygen, contact the supplier or the manufacturer.

** - 5 liters Concentrator

Here are some conditions that may require supplemental oxygen, either temporarily or long-term:



COPD



Post Surgery



Pneumonia



Sleep Apnea

Manufactured by:

Deck Mount Electronics Pvt. Ltd 260, Phase IV,
Udyog Vihar, Sector 18, Gurugram, Haryana-122015

☎ 1800 30 92 499

✉ wecare@deckmount.in

🇮🇳 Made in India

🌐 www.deckmount.in