





No Internet required


No battery needed



Free demo by experts


99.7% trace accuracy


Free ECG Interpretation

Rhythm ECG Device

Rhythm ECG device and Rhythm software supports Android. Designed for Home, Hospitals and Clinics





Lead II ECG

Lead II based Arrhythmia test for initial detection of 3 different kind of Heart Rhythms.



7 Lead ECG

Lead II based Arrhythmia test for initial detection of 7 different kind of Heart Rhythms.



Key Features

- CDSCO Certified
- One can easily record, save and share unlimited ECG reports via your smartphone
- Compatible with all android devices
- No Internet or battery is required
- Printable ECG reports on A4 size papers
- ECG trends analyzing
- Best fit for the emergency and regular purpose
- Use of chest leads for better ECG tracing

Technical Specification

Commodity Name	Rhythm ECG
Colour	White and Red
Temperature	5° C to 55° C
Humidity	15% - 80%
Model Number	Rhythm V1.0
Dimensions	64x43x12mm
Compatibility	Supports all Android Devices
Input Voltage	5V
Water Resistance	No
Warranty	1 Year
Report	PDF
Package Content	1 x ECG, 1 x Chest Lead, Electrodes, 1 Connector
Connection Type	Type C

*Term & Conditions Apply

Product Code

RHYTHM

RM101



wecare@deckmount.in



www.deckmount.in



Toll Free No.
1800 30 92 499



Deck Mount Electronics
Pvt. Ltd, 260, Udyog Vihar,
Phase 4, Gurugram 122015,
Haryana, India



Make In India