



8 Channel - Sleep Test Device

For Adults/ Paediatrics

Better User Experience | Light Weight | Easy to Wear
Real Time Waveform | Rechargeable Battery | Environmental Friendly



Friendly Interface & GUI design

Consistency | Clarity | Simplicity
Focus on the most important parameters.



Rechargeable Lithium Battery

Economic and environment friendly, supportable for 16 hours continuously monitoring.



Key Features

- Smart Appearance, only 164 g with battery
- Easy to use
Quick to learn, high reliability, low maintenance cost.
- 2.8" color screen
Real time display of all waveform data, intuitive and efficient.
- Intelligent analysis software
Automatically analyze all kinds of breathing events and accuracy is above 95%,very useful for Healthcare providers, Nursing staff and Physicians.
- Multiple sleep monitor reports, facilitate physician know patient's condition, and formulate a reasonable therapeutic schedule.
- Convenient user data management.
- Original data pages to print.

Product Code

Sleep Sense

SS100

Technical Specification

9 Parameters	Airflow, Snoring, Spo2, Chest/ Abdomen Movement, Pulse Rate, Pleth Waveform, Body Movement, Body Position.
Data Resolution	AD : 12Bit
Data Sampling Rate	50-200 Hz
Storage Data	8 GB, Over 16000 hours data
Communication	SD Card
SpO₂	70 - 100%
Pulse	30 - 240 bpm
Battery	Integrated Rechargeable Lithium Battery 4000 mAh
Continuous Recording Time	16 Hours
Screen Size	2.8 inch TFT, 320x240
Device Net Weight	214 g with battery
Gross Weight	637 g
Device Dimensions	125mmL x 76mmW x 25mmH
Report	PDF Format

*Term & Conditions Apply



wecare@deckmount.in



www.deckmount.in



Toll Free No.
1800 30 92 499



Deck Mount Electronics
Pvt. Ltd, 260, Udyog Vihar,
Phase 4, Gurugram 122015,
Haryana, India



Made In India