

Intelligent Breathing BiPAP

*Breathing Innovation,
Personalized for You!*



RESPIRARE VT 300 *ST*


A top-quality BiPAP machine delivers reliable, adjustable pressure levels with advanced comfort features and quiet operation for effective respiratory support.

BiPAP VT 300ST



Alarm Function

The alarm level meets national standards, sound alarm + graphic alarm, timely remind the cause of the alarm, for the safety of users.



Technical Specification

Modes	CPAP, S, T, ST, APCV, VAF-ST, VAF-APCV
Operating/ Monitoring Range	CPAP :4-20 cmH2O IPAP: 4-25 cmH2O EPAP: 4-23 cmH2O VT: 200-2000 ml
Sensitivity	I sens - Auto, 1-6 level, E sens - Auto, 1-6 level
I Slope	1-6 level
Ramp Time/ I Slope	0-60 mins/ 1-6 level
Insp time/ Max Insp time/ Min insp time	0.5 - 4s/ 0.5 - 4s/ 0.5 - 4s
BPM	3-40 BPM
Alarm	Power failure, Low voltage, Obstruction, SD card full, High Leak
Humidification Level	0-5 level
Display	2.8 in LCD
Dimensions (LxWxH)/ Weight	250x181x114 mm/ 1340 g
Data Management	1/7/31 days, 6 months, 1 year of detailed data.
Power Adapter	Input Range : 100-240V, 50/60 Hz, 1.8 A max (typical power consumption)
Operating Temperature/ Humidity/ Altitude	+5 °c to +35 °c/ 10-95% relative humidity, non-condensing/sea level to 2,591 m, air pressure range 700 to 1060 hpA
Storage and Transport Temperature/Humidity	-20 °c to +60 °c / 5-95% relative humidity, non-condensing
Supplement Oxygen	Maximum Flow: 4L/min (CPAP), 15 L in S, T, ST, APCV, VAF-ST, VAF-APCV
Standard Air Filter	Cotton Filter
Tubing/ Air Outlet	Standard (19 mm)/ 22 mm Air Outlet (STD)/ Heating tube (Optional)
Sound Pressure Level	<30 dBA with uncertainty of 2 dBA
Warranty	2 Years (Repair)
Conforms to Standards	IEC 60601-1-1 Class 2 (double insulated, IEC 60601-1-2 Type BF, Ingress Protection IEC 80601-2-12 IP22) CDSO

*Term & Conditions Apply

Key Features

Auto On/Off

This feature will turn the device **On** for you have the mask on without the need to press any buttons and will also turn the device **Off** once you take your mask Off.

Integrated Humidification

Heated Humidification makes therapy more comfortable and can reduce side effects like a dry nose and throat.

Mobility, Portability

Easy to carry the device from one location to another.

Product Code

Respirare BiPAP with Humidifier **VT 300 ST**

Note: Technical specification may change without notice. Always refer to your device user manual.

Deck Mount Electronics Pvt. Ltd

260, Phase IV, Udyog Vihar, Sector 18, Gurugram, Haryana 122015

1800 30 92 499

Made in India

wecare@deckmount.in

www.deckmount.in